

3D CELL CULTURE SCAFFOLDS

GROWTH AND INTERACTIONS IN ALL THREE DIMENSIONS

Product concept

3D cell culture scaffolds designed as two-layer constructs:

- The top layer features a unique matrix of fibres and pores fabricated by the proprietary 3D fibre printing method, and serves as a matrix hosing cells.

- The bottom layer features a dense network of nanofibres forming a cell impenetrable membrane, preventing cell migration to the well floor.

Such design provides a favourable environment for cell attachment and proliferation. The round fibre and pore morphology follows a random structure of natural extracellular matrix (as compared to straight/aligned fibre and filament structures). Large pores allow for an efficient cell distribution throughout the entire scaffold height.

Unique and Fine-tuned Scaffold Morphology



3D-PCL/PT-08/035-NF Fine fibres and small pores 3D-PCL/PT-22/125-NF Medium fibres and pores 3D-PCL/PT-25/155-N Coarse fibres and large pores Nanofibrous layer for cell support

Bious[®] 3D cell culture scaffolds

Specifications:

Construction: 3D, 2-layer (nanofibre bottom) Basis polymer: $Poly(\epsilon$ -caprolactone), plasma-treated

SKU: **3D-PCL/PT-08/035-NF** Fibre diameter (mean) 8 μm Pore diameter (mean) 35 μm

SKU: **3D-PCL/PT-22/125-NF** Fibre diameter (mean) 22 μ m Pore diameter (mean) 125 μ m

SKU: **3D-PCL/PT-25/155-NF** Fibre diameter (mean) 25 μm Pore diameter (mean) 155 μm





Confocal microscopy images of MDA-MB-231 cells in scaffold

Features

- Tunable morphology
- · Cut-to-size
- Easy imaging
- · UV-sterilized and ready to use
- Compatible with most current 2D assays

Applications

- · Drug screening
- Environmental cytotoxicity
- Tissue engineering
- Organ-on-chip models

Modifications

- Wetting properties
- Cold plasma surface activation
- Attachment of functional molecules (growth factors)



3D CELL CULTURE SCAFFOLDS

GROWTH AND INTERACTIONS IN ALL THREE DIMENSIONS

Selectable options

Scaffold diameter21, 15, 5 mm (fits 12, 24, 96 well plates, possible to use w or w/o support holders)Scaffold thickness $200 \pm 20 \ \mu m$ Fibre diameter $10 - 200 \ \mu m$ Pore diameter $15 - 500 \ \mu m$ Morphologyround interconnected poresWater contact angle $40^{\circ} - 70^{\circ}$

Method work flow



- 1. Unpack test-plates in a laminar flow box and inspect visually.
- 2. Fill wells with 70 % ethanol for 20 min and then aspirate. Wash with sterile PBS (three-times).
- 3. Pre-wet wells with culture media: add culture media, incubate for at least 30 min and aspirate the medium.
- 4. Prepare cell suspension in a required density and pipette the appropriate volume of cell suspension (refer to 3D cell seeding protocol).
- 5. Carefully place cell suspension droplet onto the centre of the well surface.
- 6. Allow droplet to distribute through the entire scaffold.
- 7. Refill the wells with culture media and change the media every 3 4 days.
- 8. Perform selected analytical procedures.

Common protocols for 2D cultures may be adapted to 3D cultures. Scaffolds are compatible with both assays.

Recent application: https://doi.org/10.1016/j.bej.2022.108531



Supported by the Lithuanian Innovation Agency under the project "Hybrid 3D scaffolds for in vitro tissue engineering research", No. 09.3.3-LMT-K-712-19-0181, in cooperation with Lithuanian University of Health Sciences.



Bious Labs Street address: Draugystes g. 15D LT-51227 Kaunas Lithuania Web: life.biouslabs.com Email: sales@biouslabs.com Phone/WhatsApp: +370 675 24358 Bious Technology & Life